

# Sand Baths



GHP-9S  
300X300mm

## SHP-4X2D/HP-400S, DIGITAL SAND BATHS

- ▶ Robust construction
- ▶ Eloxidised aluminium heating surface
- ▶ Good uniformity over the entire heating surface, including the corners and edges
- ▶ Temp. can be displayed in °C or °F.
- ▶ Temp. stability:  $\pm 1^\circ$
- ▶ PID temperature controller
- ▶ Digital set-point & read-out of actual temperature
- ▶ 220Volt/50Hz 3400watt.

MODEL	GHP-9S	SHP-4X2D	HP-400S	HP-600S
Working range	up to 400 °C			
Stability	$\pm 1^\circ$			
Setting & Readout	Digital			
Heater	1600W	3000W	3400W	4500W
Top-Plate(mm)	300x300mm	250x500mm	400x400mm	400x600mm
Dimensions WxHxD (mm)	300x220x300	500x220x250	400x270x400	400x270x600
Power source	220V 50Hz			



HP-400S  
400X400mm



HP-600S  
400X600mm